



NOUVEAU CONTOUR®
MASTERS IN MICROPIGMENTATION

Needles



Safety Cartridge Needles

Nouveau Contour is extremely concerned with the threat of fluid cross-contamination to the hand piece. The Nouveau Contour patented, Safety needle system prevents from fluid cross contamination. Verifiably, traditional machines are subject to airborne and fluid cross-contamination.

Unique Needle Cartridge System

Tests have proven that the risk of cross contamination with traditional machines due to possible backflow of liquids and circulation of contaminated air within the hand piece is prevalent.

For the safety of both the technician and the client, the Nouveau Contour Needle Cartridge System works unlike anything on the market today. The needle and the needle "tube" are integrated in one disposable cartridge system. Internal diaphragms halt and prevent fluid and airborne contaminants from entering the hand piece, which is therefore completely independent of the disposable cartridge and void of any possibility of cross-contamination. When the cartridge is detached from the hand piece the needle automatically retracts back into the cartridge, preventing an accidental "stick" to the technician, or client. The safety needles fit on the Nouveau Contour safety hand piece. The system is therefore void of any possibility of cross contamination. Available in 21 different needle formations.

Safety first!

Nouveau Contour realizes that it is imperative that the disposable tube/needle cartridges be sterilized. Each cartridge is packed separately in a class 100.000 clean room and sterilized using EO gas to prevent disease transmission. The entire cartridge, which contains the needle, can then be disposed of properly. Another safety feature is that when the cartridge is detached from the hand piece the needle automatically retracts back into the cartridge, preventing an accidental "stick" to the technician, or client.

Safety Needles



Round-tip Needles

1-point Nano NT needle 0.12MM

Ultra fine, high quality needle. The steel tube at the tip gives less vibration on the needle. The innovative tunnel provides more stability which results in extraordinary pigment deposit and ultra sharp hairstrokes. Only for experienced technicians.



1-Point Nano 2 Needle 0.12 MM

Create ultra-sharp and crispy lines with our thinnest nano needle. Generates smaller size perforations. Always use 2 tones darker than usual for optimal results. For experienced technicians only.



1-Point Nano 1 Needle 0.15MM

Our second thinnest needle, this penetrates the skin shaft at 0.15MM. Always use 2 tones darker than usual for optimal results. For experienced technicians only.



1-Point Micro Needle 0.18MM

Great for creating fine hair strokes. Be careful not to go too deep into the skin as there is less resistance of the skin when using the micro needle.



1-Point Liner Needle 0.25MM

Suitable for fine hair strokes and detailed lines, such as lip corners. Remain watchful against signs of scarring and do not inject the pigment too deeply to avoid too dark a color.



3-Point Micro Needle 0.18MM

The 3-point micro needle is thinner than the regular 3-point needle. Suitable for ultra-fine lines, such as lower eyeliner and lip corners. Do not inject the pigment too deeply to avoid too dark a color.



3-Point Outline Needle 0.25MM

A 3-point needle whereby one of the points protrudes in the direction of the supply opening. Suitable for making a distinct contour to define a high eyeliner, semi hair strokes, and the outer corners of the lips.



3-Point Liner Needle 0.30MM

Suitable for drawing lines in eyebrows, fine lip lines, and eyeliner with lash-broadening effect.



5-Point Round Needle 0.30MM

Suitable for soft eyebrows, wide eyeliners and lip liners.



5-Point Shader Needle 0.30MM

The 5 shader is commonly used to create soft shaded brows, fill in lips and to create a smudged look for eyeliner. It is also a good choice for camouflage work and to fill in the areola.



7-Point Round Needle 0.30MM

Has an open, round needle form and is suitable for very wide upper eyeliners and the partial coloring in of the lips.



Flat-tip Needles

4-Point Flat Needle 0.35MM

Four needles in a row. Suitable for the hair-stroke technique for eyebrows, by applying the 'scoop' technique. The result is a very natural look of eyebrow hairs. (Master Training Course Hair stroke)



5-Point Magnum Needle 0.30MM

Flat form comprising 5 needles in 2 rows. Suitable for smokey eyeliners, wide and high eyeliners, and the coloring in of lips using the circling technique. (Master Training Course Advanced Eyeliner) (Master Training Course Advanced Lip)



9-Point Magnum Needle 0.30MM

Flat form comprising 9 needles in 2 rows. Suitable for a subtle coloring in of large areas. (Master Training Course Areola) (Master Training Course Advanced Lip)



Digital Blading Needles

Digital Blading 3-Point Sloped-Tip Needle 0.40MM

Flat form comprising, 3 nano needles in a row. For experienced technicians only. Designed especially for digital blading. Creates razor sharp hair strokes and is suitable for all skin types. Suitable for 'digital blading' in a somewhat thinner eyebrow.



Digital Blading 5-Point Sloped-Tip Needle 0.40MM

Flat form comprising, 5 nano needles in a row. For experienced technicians only. Designed especially for digital blading. Creates razor sharp hair strokes and is suitable for all skin types. Suitable for 'digital blading' in a somewhat bolder eyebrow.



Sloped-tip Needles

3-Point Sloped-Tip Needle 0.40MM

Flat form comprising 3 needles in a row. Suitable for relatively high eyeliners, a wide lip liner, or the partial coloring in of lips. (Experience desired when working with this needle)



5-Point Sloped-Tip Needle 0.40MM

Flat form comprising 5 needles in a row. Suitable for a very high, advanced eyeliner, and full lips without lip liner. (Master Training Course Advanced Eyeliner) (Master Training Course Advanced Lip)



Power-tip Needles

3-Point Power Needle 0.30MM

Extra sharp 3-point needle for more precise work than a regular 3-point needle will allow. Suitable for the shading of eyebrows and the drawing of fine eyeliners and lip liner.



5-Point Power Needle 0.30MM

Extra sharp 5-point needle for more precise work than a regular 5-point needle will allow. Suitable for the shading of eyebrows and the drawing of somewhat wider eyeliners.



7-Point Power Needle 0.30MM

Extra sharp 7-point needle for more precise work than a regular 7-point needle will allow. Suitable for the subtle shading of eyebrows and the drawing of wide upper eyeliners, a wide lip liner, and a partial shading of lips.



Needle formations

1-point nano NT	•	3-point micro	••	5-point shader	••••
1-point nano 0,20 mm	•	3-point power	•••	5-point power	•••••
1-point nano 0,25 mm	•	3-point sloped		DB 5-point shader	
1-point micro	•	DB 3-point sloped		5-point sloped	
1-point liner	●	4-point flat	••••	7-point round	••••••
3-point liner	•••	5-point round	••••	7-point power	••••••
3-point outline	•••	5-point magnum	••••	9-point magnum	••••••••



Contact us!

Nouveau Contour Headquarters

Burgemeesterlaan 2 | 6002 EG Weert | The Netherlands

P. +31 88 - 51 55 810 | **E.** export@nouveaucontour.com | www.nouveaucontour.com